

RediHeal - Which Version is Correct for You?

	RediHeal Wound (Dry Fibers)	RediHeal Wound (Ointment)	RediHeal Dental
Soft tissue promoter / wound healing	X	X	X
Healthy bone promoter			X
Wound Types:			
Heavily effusive wounds	X		
Lightly effusive wounds (bandageable)	X	X	
Lightly effusive wounds (non-bandageable)		X	
Dry wounds	X (with wet/dry)	X	
Degloving injuries	X		
Mass cell tumor removals (non-healing)	X	X	
Bite wounds (non-puncture type)	X		
Minor cuts, scrapes and scratches		X	
Burns	X	X	
Non-healing / chronic wounds	X	X	
Dirty wounds / Abscesses	X		
Dental extraction sites			X
Periodontal pockets			X
Dental abscesses	X		X
Oronasal fistulas	X		X
Suitable for:			
Equine	X	X	X
Companion animals	X	X	X
Avians / exotics	X	X	X
Reapplication	app. 3-7 days	Daily (minimum)	Typically none required

NOTE: This material is not recommended for use as a hemostat; it may be used in conjunction with a hemostat for significant bleeding. It will provide very light hemostatic benefits.